Additional information
Zusatzinformation
Informations supplémentaires
Información adicional

Additional information

Model PG43SA-C, autoclavable version

EN

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Prior to starting any work, read the operating instructions!
 Keep for later use!

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
Supplementary documentation:

- ▶ This additional information applies in conjunction with the operating instructions “Pressure gauges for sanitary applications, models PG43SA-S, PG43SA-C” (article number 11532319).

This additional information supplements the respective chapters of the above-mentioned document.

2. Design and function

2.3 Vent valve

By opening or closing the vent valve, the case will be vented to or sealed from the atmosphere. Only ever operate the vent valve via the spanner flats and to the prescribed torque using an appropriate tool. By opening the vent valve, the pressure difference between the atmosphere and the sealed case will be eliminated. Maintaining the ingress protection is only ensured with a closed vent valve.

5. Commissioning, operation

EN

5.5 Commissioning and operation

For pressure compensation, open the vent valve for a short time and then close it again directly, see chapter 2.3.

The vent valve remains closed during operation.

Only through a brief venting of the case before each pressure measurement can the pressure be read correctly.

5.7 Autoclaving

- Ensure that, before the start of the autoclaving process, the vent valve has been closed.
- Carry out the autoclaving process. For performance limits see chapter 9 “Specifications”.
- Allow the instrument to cool to ambient temperature.
- By briefly opening the vent valve, the pressure difference between the atmosphere and the sealed case will be eliminated.

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All rights reserved. / Alle Rechte vorbehalten.
WIKA® is a registered trademark in various countries.